



MILESTONE
H E L P I N G
P A T I E N T S

EasyCUT

Static Remover for Microtomes

REMOVE STATIC CHARGES FROM ANY PARAFFIN BLOCK PROCESSED

During paraffin block cutting, electrostatic charges may often occur. These charges make ribboning difficult, as the sections stick to the paraffin block itself, to forceps and to user's fingers. Transfer to the water bath is almost impossible and section ribbons may break apart.



Static charges may result in:

- **Delays** in slide preparation
- **Longer** microtome cutting sessions
- **Stressful** workload for the operator
- **Unconfident** user procedure

EasyCUT Benefits and Values:

EasyCUT neutralizes electrostatic charges and thus helps technicians in obtaining a perfect, non-sticky, wrinkle-free ribbon, without shattering effect, which is easy to transfer to the water bath.

SMOOTHER CUTS

For more efficient slide preparation

FASTER BLOCK SECTIONING

Improves efficiency and turn around time. Its built-in led lamp provides ergonomic working conditions in low light.

WRINKLE-FREE SECTIONS

The section ribbons are perfectly stretched with neither shattering nor scattering. They do not break up in the water bath.

EasyCUT FITS ON ANY ROTARY MICROTOME

Its universal structure makes it fully compatible with any laboratory microtome set up.

NO EXTRA BENCH SPACE NEEDED

EasyCUT only requires the footprint of microtome and it does not impair routine procedures.

24/7 OPERATIVITY

Working non-stop. Soundless.

INDUSTRIAL GRADE APPLIANCE

It's made to last. EasyCUT is resistant to laboratories' harsh conditions for long lasting durability, with no need for maintenance.

EasyCUT Main Features:

- Air fan with **4 speeds**
- Built in **LED** lamp
- Air **filter**
- **Soundless** operations
- EasyCUT adapts to **any** rotary microtome on the market, both manual and automatic
- It does not occupy extra **bench space**
- **Maintenance free**

Headquarters

MILESTONE Srl
Via Fatebenefratelli, 1/5
24010 Sorisole (BG) - Italy
Tel: +39 035 4128264
Fax: +39 035 575498
email: marketing@milestonemed.com

US Headquarters

MILESTONE MEDICAL TECHNOLOGIES, INC.
6475 Technology Avenue, Suite F
Kalamazoo, MI 49009 - USA
Tel: 269-488-4950 - Toll-free: 866-995-5300
Fax: 269-488-4949
email: info@milestonemed.com