Seeing Beyond Limits




HLA Expertise, Simplified.

The MATCH IT![®] Antibody software suite was designed specifically as an accessory to assist in the evaluation of test results from the LIFECODES Antibody portfolio.



MATCHIT!

Introducing MATCH IT! Antibody, a software designed specifically for simplifying the complexities of HLA.




Streamlined Data Analysis


- Built in CREG and Eplet Analysis tools
- Advanced highlight and overlay tools to help determine antibody specificity
- Multiple Data Graphing features
- Automatic import of Luminex data
- Options for server based or local database installation
- Remote support provided by Immucor support teams
- Multiple custom utilities available to manage import of data into various LIS

ADAPTIVE

CUSTOMIZABLE page layouts, exports and reports fit everyone's needs!



Flexible display options that include both a **BEAD RANKING GRID AND HISTOGRAM**. Both the bead rank and histogram can be customized.









EFFECTIVE

EPITOPE ANALYSIS features that predict most probable eplet. Bead Rank report also displays **RELATIVE ANTIGEN DENSITY** for each antigen as a tool to explain reactivity.




Multiple analysis preferences such as **MFI CUT OFF** to adjust cutoffs based on the knowledge of an experienced HLA technician.



POWERFUL

PATIENT VIEW displays all batches for a given sample. They can be **GRAPHICALLY COMPARED** from one batch to another, and **ALL ANTIBODY HISTORY** viewed.



COLOR CODED ANTIGENS or groups of antigens based on serological cross-reactive groups (CREGs) for easier viewing.








VALUABLE


DIRECT C3d AND LSA IgG COMPARISON

with graphing tools and export options.



OVERLAY AND GRAPHING OPTIONS

to review data from same time points and/or for comparisons between the ID and LSA.





To learn more or to set up a demo with one of our Transplant Business Managers, **855.IMMUCOR (855.466.8267)** or visit **www.immucor.com**.

International Contact: imarketing@immucor.com



www.immucor.com

